

*Summon help instantly
in a difficult or
emergency situation:
whether or not you are
conscious or able.*

*Peace of mind for anyone – and those who
love or care for them.*

*"This tiny device is one of the
most amazing inventions so far
this millennium. It's so simple
to use yet its capabilities are
easily life-saving. Every potentially
vulnerable person should have
one – and its additional **tap2tag®**
capability is essential for anyone with a
diagnosed medical condition."*

Dr R Hopewell, GP



Pebbell+ powered by **tap2tag®** enables instant communication between a user and their emergency contacts, even alerting them to a fall or wandering. Find out how...

Pebbell+ powered by **tap2tag®** is a mobile pendant alarm – a two-way communication device connecting a potentially vulnerable person with someone who can help in an emergency or difficult situation.

It's the only mobile pendant alarm in the world powered by **tap2tag®**, meaning it's also equipped with NFC (near-field communication) technology. This gives first-responders and paramedics access to vital personal medical information such as known conditions, prescribed drugs and life-saving interventions.



How does the device work?

It's a small, discreet aid, about the size of an egg. A typical user wears the device around their neck or attached to a belt, although it can be carried in a pocket.

It's very simple to use, yet has incredibly advanced features and capabilities.

Fully waterproof : stay safe, even in the shower!



PEBBELL+

POWERED BY

tap2tag
Medical Alerts®



Two-way communication: **or one-way in an** **emergency**

At its most simple, it's a tiny mobile phone with a SIM card, microphone and speaker, which users and their carers or loved ones can use to communicate.

In a difficult situation or an emergency, however, it becomes a potentially life-saving device. If the wearer is in difficulty, they simply press a button to summon help. If they have a fall, the device detects the sudden drop, through an in-built accelerometer, and alerts emergency contacts. You can pre-programme up to three emergency contacts into

the device. In the event of activation or a fall, the device sends a text message (SMS) to each contact giving:

- details of the device owner (eg. 'Mum', 'Dad', 'Joan')
- date and time of activation
- a link to Google Maps™ identifying the exact location of the user
- 'fall alert' if appropriate
- current battery life of the device

The device also telephones the first-named contact. If the call is not answered it calls other contacts. If the patient is conscious they can speak to the emergency contact.

Life saving information... *on tap.*

Wander alert!

Pebbell+ powered by tap2tag® is also equipped with the latest 'geo-fencing' software. This is a powerful feature, extremely helpful in preventing vulnerable people from wandering beyond where they are safe, such as their home and garden or the boundaries of a care home.

You can specify the area outside of which the alert activates, for example beyond 50 metres of their home. If the 'wander' alert is activated the text messages also give the exact location of the wearer. At any time, a carer or loved one can text a code to the device and it will return information on their whereabouts.

Life-saving information on tap: the only pendant alarm embedded with tap2tag®

Your device also includes all the benefits of tap2tag® Medical. You can upload vital medical information to your Pebbell+ which becomes immediately

accessible to first-responders, family or passers-by in the case of an emergency.

You choose the data you wish to share. Typical users include:

- their name
- known medical conditions
- any and all medications they need
- any message, for example

"Call 999. I have a heart condition. Spray is in my bag. Spray once in my mouth then wait for paramedics."

A user might also choose to include information on organ donation, their date of birth, blood group and a photograph.

No data is stored on the device. You upload the information online on a very secure, password-protected portal, and you can update it at any time. It's a great way of keeping track of medications!

How do first-responders access the information?

Every **tap2tag®** device is embedded with ultra-modern 'NFC' technology (near-field communication). It is used through a mobile phone. In an emergency, for example if the wearer collapses, after calling the emergency services anyone with an NFC-enabled mobile phone can 'tap' the device. With a couple of clicks they can gain instant access to the vital information. Even if the first-responder hasn't activated the **tap2tag®** Medical information, the paramedics can do so on arrival, gaining potentially life-saving information about the patient.

All modern phones have NFC, but in the unlikely event that the first on the scene doesn't have NFC they can access the information using any internet-enabled device, including their mobile phone.

In the meantime, of course, the three contacts pre-programmed into the device have received text messages and phone calls alerting them to the incident.

"This little device is a real gem. I can't begin to understand how it works, but it does. My daughter and son-in-law bought one for me two months ago and it's lifted a great weight off my mind – and theirs!

I think it provides personal security in many different ways for any potentially vulnerable person and relaxes those who care for them. You don't have to be very old or very ill to benefit enormously."

*Mrs Wakley
Bristol*



You can link your Pebbell+ powered by tap2tag® to a call centre careline such as Centra Pulse or a local authority service. It can be your first-named contact if you prefer.

How much will I pay for a Pebbell+ powered by tap2tag®?

Perhaps not as much as you might think, for a device which gives you and your loved ones peace of mind and which may potentially save your life!

There is a one-off charge of just £168 for the device, including :

- Charger
- Stand
- Highly durable cord attachment
- £8 credit for texts and phone calls.

After that there's just a small monthly fee for texts and phone calls, from £3 per month plus calls and text messages.

How to buy your device

You can buy a Pebbell+ powered by **tap2tag®** device from our secure website at www.tap2tag.me. The account set-up process is very simple.

PEBBELL+

POWERED BY **tap2tag**
Medical Alerts®

*Available in four modern colours:
blue, pink, lime & black.*



Visit www.tap2tag.me for more information.

Pebbell+ powered by **tap2tag®** and **tap2tag Medical** were developed and are promoted by tap2tag Limited. Registered offices at 41 High Street, Kingswood, Bristol BS15 4AA, United Kingdom